

Spor-Klenz Concentrate Cold Sterilant is a fast-acting liquid sporicide specifically formulated for sterilization and disinfection of environmental, hard, non-porous surfaces. Spor-Klenz Concentrate Cold Sterilant contains the proven blend of peracetic acid, hydrogen peroxide and acetic acid to provide fast, effective microbial control of challenging spores, such as *Bacillus subtilis* and *Clostridium sporogenes*. This formulation can be diluted with water, offering operational flexibility for use as a disinfectant or sanitizer. This product has an extensive documentation package to meet your validation objectives and is supported by the STERIS Technical Services team for individualized application and validation assistance.



## Spor-Klenz™ Concentrate Cold Sterilant

Fast-acting sporicide in a versatile format

### Features

Concentrated, easy to dilute formulation

Broad spectrum efficacy against bacteria, fungi, viruses and bacterial spores using AOAC protocols

Excellent material compatibility

Low toxicity profile

0.2 micron filtered and each lot USP sterility tested

Manufactured in accordance with cGMP with complete lot traceability

### Benefits

Can be used as a sporicide, disinfectant or sanitizer

Reduces packaging waste and required storage space

Confirms efficacy against many common microbes

Suitable for most surfaces


Supports operator safety

Facilitates introduction into critical areas

Supports quality assurance and quality control goals

### Applications

- Surface disinfection and microbial control of:
  - » Floors
  - » Walls
  - » Ceilings
  - » Countertops
  - » Biosafety Cabinets
  - » Carts and Other Equipment
- Product may be applied by:
  - » Mopping
  - » Spraying
  - » Wiping



<sup>1</sup>Products that fall into the categories at the bottom of the pyramid are most frequently used and are generally not sporicidal. Progression up the pyramid indicates stronger performance overall and a broader spectrum of claims.

In maintaining microbial control of your facility, sporicides perform a critical function, destroying the toughest environmental contaminants. They should be used based on the frequency of spore results in your environmental data and in rotation with one or more surface disinfectants. Spor-Klenz Concentrate Cold Sterilant offers fast, broad spectrum efficacy with superior material compatibility. Our Technical Services team can help customize your contamination control program to take advantage of this excellent product.

## Ordering Information

### Spor-Klenz Concentrate Cold Sterilant

Filtered for consistent quality

Description		United States	Canada
32 oz (946 mL) bottle	4 per case	652026	652026CA
1 gallon (3.8L) bottle	4 per case	652008	652008CA
40 gallons (151.4L)	55 gallon drum	6520T7	6520T7CA